

## 2-4 Micro-plate Centrifuge

Marks Scientific



2-4 micro porous plate centrifuge is our company specializes in 96-hole design of the instantaneous centrifuge to facilitate the liquid separation from the wall. Most micro plate centrifuges on the market are bulky and occupy a large space of laboratories, this micro plate centrifuge is unique designed and small, only 23x20cm. The micro plate is vertically loaded into the rotor from the top slot of the centrifuge, and the liquid is kept at the bottom of the micro plate due to the surface tension. Therefore, it will not leak at all

### **Key Features:**

- ⦿Rapid centrifuge hanging wall drop;
- ⦿Used for PCR experiment effectively;
- ⦿Multiple accessories included PCR micro plates
- ⦿Less than 1/4 volume of standard micro plate centrifuges in the market

**Technical parameters:**

<b>Model</b>	<b>Micro-plate Centrifuge 2-4</b>
RPM (rpm)	2500rpm
RCF	600×g
Weight	2.1KG
Capacity	96 Hole PCR Plate
Voltage	AC120V/220V,50/60HZ
Dimensions	190x220x186mm