

L3-5KM / L4-5KM Beauty centrifuge

Platelet Rich Plasma Centrifuge , Self-fat transplantation centrifuges:



Features:

- ⦿ Medium size with high compatibility.
- ⦿ Microcomputer control with LED display, driven by AC frequency conversion motor.
- ⦿ Special design with Shock absorber ,auto-balance function.
- ⦿ Electronic locks with full steel inner cavity protective cover.
- ⦿ RCF and speed separately and switching free, automatically calculate RCF.
- ⦿ Automatic calculates and displays RCF at the same time.
- ⦿ 10 level speed control.
- ⦿ One key short-time centrifuge.

Technical Parameter:

Model	L3-5KM	L4-5KM
Max RPM	4000rpm	4000rpm
Max RCF	2750×g	2750×g
Capacity	10ml、20ml	10ml、20ml、50ml
Timer	1min~99min	1min~99min
Revolutions/min	±10r/min	±10r/min
Voltage	AC 220±22V 50Hz 10A	AC 220±22V 50Hz 10A
Power	450W	750W
Noise Level	≤ 60dB (A)	≤ 60dB (A)
Chamber Diameter	Φ380mm	Φ420mm
Outer Dimensions	525×430×360 (mm)	560×460×400 (mm)
Packing Dimensions	600×500×440(mm)	640×530×470(mm)
Net Weight	40kg	50kg
Gross Weight	48kg	55kg

PRP injection transplantation centrifuge:

In order to fully discharge PRP, PRP injection and beauty centrifuge has done a lot of research on rotor, rotation speed, centrifugal force and lifting speed. It has greatly improved the effective extraction rate of PRP and shortened the operation time. It can be equipped with a special PRP kit imported from South Korea, which greatly improves the effective extraction efficiency of PRP, making the whole treatment process easier, faster and more effective.

Autologous fat transplantation Centrifuge:

Autologous fat transplantation beauty centrifuge is mainly used in plastic surgery, dermatology and other fields. Using 10ml, 20ml, 50ml conventional syringes, rapid separation and purification of fat, no fever, no damage to fat cells, greatly ensuring the survival rate of fat transplantation, reducing the complicated process of transplantation, and ensuring effective fat in fat transplantation And fat survival rate .