

L4-4KR Floor Low Speed Refrigerated Centrifuge



Key Features:

- ⦿ Microcomputer control with LED display, driven by AC frequency conversion motor.
- ⦿ RCF and speed separately and switching free, automatically calculate RCF.
- ⦿ Up to 10 level accelerate and decelerate control.
- ⦿ Over speed, over heat, unbalance protections with auto lid lock.
- ⦿ Electronic locks with full steel inner cavity protective cover.
- ⦿ Emergency lid safety lock.
- ⦿ Rotor parameter locking function.
- ⦿ Imported Eco-Friendly refrigerating system.

Applications:

L4-4KR Floor Low Speed Refrigerated Centrifuge is widely used in hospitals, universities and colleges at all levels due to its high performance and reliable quality.

Technical Parameter:

Model	L4-4KR
Max RPM	5000rpm
Max RCF	4390×g
Max Capacity	4×750ml
Timer	1min~99min59s
Revolutions/min	±10r/min
Temp Range	-20℃~40℃
Temp Accuracy	±1.0℃
Voltage	AC 220±22V 50Hz 15A
Power	1300W
Noise	≤ 60dB (A)
Chamber Diameter	Φ400mm
Dimensions (L x W x H)	680×580×850 (mm)
Packaging dimensions	780×680×950(mm)
Net Weight	100kg