

---

## L6-6K Floor Low Speed Large Capacity Centrifuge



### Key Features:

- ⦿ RPM stepped in 100rpm with accuracy  $\pm 10$ rpm.
- ⦿ Acceleration/deceleration time customized.
- ⦿ Auto detection of imbalance, over speed, over heat and door lock protection.
- ⦿ Self-correction balancing swivel (Optional), up to 3 \* 50g weight imbalanced.
- ⦿ Full stainless steel body ,more stable and quite during operation ,noise level  $\leq 60$ dB(A).
- ⦿ Brushless AC frequency conversion motor , maintenance free.
- ⦿ High capacity up to 6 \* 1000ml swing out rotor ,especially for blood bank.

### Applications

L6-6K has large capacity, stable and durable, with continuously optimized electrical control, specially designed for customers with large capacity requirements, such as blood stations, pharmaceuticals, bioengineering and etc.

**Technical Parameter :**

Model		L6-6K		
Max RPM		6500rpm		
Max RCF		6680×g		
Max Capacity		6×1000ml		
Timer		1min~23h59min		
Revolutions/min		±10r/min		
Voltage		AC 220±22V 50Hz 30A		
Power		1500W		
Noise Level		≤ 60dB (A)		
Chamber Diameter		Φ600mm		
Outer Dimensions		840×730×1240 (mm)		
Packing Dimensions		940×830×1340(mm)		
Net weight		190kg		
Gross weight		220kg		
Rotor	NO.	Capacity	Max Speed ( r/min)	Max RCF(×g)
Angle Rotor	1	6×250ml	6500rpm	6610×g
	2	6×300ml	6500rpm	6610×g
	3	6×500ml	6000rpm	6680×g
Swing-out Rotor	4	6×1000ml Round Type	4200rpm	5100×g
		6×1000ml Square Type	4200rpm	5100×g