

MiniStar 4k 7k 10k Table Mini Portable centrifuge



Key Features:

- ◎ Smooth and portable design with low noise level.
- ◎ Safety design with flip switch function, automatically stop while open; Cover made of high strength composite material.
- ◎ Aviation level micro DC brushless motor
- ◎ Most advanced broadband power technology designed, input voltage 85V ~ 245V
- ◎ Machine provide two rotor and two types of tubes

Technical Parameter:

Model	MiniSt 4K	MiniSt 7K	MiniSt 10K
Max Speed	4000rpm	7200rpm	10000rpm
Max RCF	2000xg	2350xg	5000xg
Noise	≤50db	≤52db	≤54db
Rotor	8×1.5/2.0ml	0.5/0.2mltube	2×8×0.2ml
Voltage	85~245V/50-60HZ		
Packing	150×180×120mm		
Weight	1.5kg		