



## 100g Universal Powder Machine Grinder MRDFT-100 T



### ***Specification: Model: MRDFT-100 (Dry Grinder)***

Voltage: 220V 50HZ. Motor power: 400 W

Applicable materials: universal

Material: Stainless steel #304

Capacity: 10-100 g

Production capacity: 1.2-1.5 (kg/h)

Speed: 25000 rpm

Fineness: 50-200mesh

Dimensions: 14\*14\*30 cm

Weight: 2.8 (kg)

Work half an hour, five-minute break

Degree of grinding: ultrafine mill

The max hardness of material: ore, pearl, etc.

**Contact: Marks Scientific Solutions Inc.,  
4021 McGinnis Ferry Road, Suite 1814A  
Suwanee GA 30024  
Phone: +1 770 814 6427  
Email: [info@marksscientific.com](mailto:info@marksscientific.com)**

**[www.marksscientific.com](http://www.marksscientific.com)**